




SECURITY GUARANTEED


100%



Anti-wrinkle removal
5D-HIFU machine



Function and Principle




Wrinkle removal



Sagging skin removal

applications




Skin rejuvenation




3D - V shape face rebuild

1. Instead of the traditional surgical face lift with a non-invasive technique;
2. With a special transmitter transmission from the skin surface conduction of high-energy focused ultrasound to the SMAS layer;
3. Proprietary high-speed transmission technology of 0.1-0.25;
4. Up to 15000 lines of a layered high-precision emitter transmitter is three times that of similar instruments, greatly improving the cost effectiveness;
5. Exclusive development of stable output 4-10MHZ broadband technology that fits well with layered skin treatment technology.
6. Quick & short treatment time: 30 MINS one face treatment.
7. 7 transducers for skin lifting and body slimming.



The treatment principle of Sincohifu System is High Intensity Focused Ultrasound(HIFU) which can be focused on a single point so as to produce high energy that acts on corium layer, SMAS layer to stimulate proliferation and recombination of collagen, then the outline will be tightened and wrinkles removed.





Handpieces Details

Face Lift Wrinkle Removal Body Slimming

- 1.5mm: to penetrate the skin to a depth of 1.5mm, is responsible for activating the skin epidermis in thinner tissue.
- 3.0mm: to penetrate the skin to a depth of 3.0mm, is responsible for activating the skin's dermal layer of collagen, effectively enhance the effect of the consolidation of the outline, but also improve large pores and reduce the appearance of wrinkles, targeted to chin, chest.
- 4.5mm: to penetrate the skin to a depth of 4.5mm, direct subcutaneous SMAS layer, targeted to thicker skin.




1 Professional HIFU Handle



- **Face Lifting**
5D hifu works for face lifting, skin tightening, wrinkles pores, V-line forming, jowl line, cheek, skin tone improvement, and laxity improvement.
- **Body Slimming**
Whole body contouring of chest, arm, abdomen, flank, thigh, lower leg, and back.


2 Gold Rf Microneedling



- New treatments double effect.
- No scab, colorless, sink, perfect effect.
- The curative effect is remarkable and the effect is maintained for 3-6 years.



3 Vmax Hifu Handle



- Eye lift & Wrinkle removal
- Smaller areas that the big 4Dhifu handle cannot do it.

Treatment Scope:

1. Remove wrinkles on around forehead, eyes, mouth, etc.
2. Lift and tighten both cheeks skin.
3. Improve jaw line, reducing "marionette lines"
4. Tighten the skin tissue on forehead, lifting the eyebrows lines.
5. Improve skin complexion, making the skin delicate and bright
6. Match with injecting beauty like Hyaluronic acid, collagen, to resolve more aging problem.
7. Remove neck wrinkles, protecting neck aging.

4 Lipo Hifu Handle



- **Fat reduction**
Equipment with 2 cartridges, reach to different fat depth.

Ultrashap+ ultrasonic focusing de-fating cosmetic instrument integrates the Ultrashap technology. Its probe can automatically scan the skin, identify lipocyte, and precisely shoot energy points into the subcutaneous fat, accumulation through non-invasive focusing ultrasound wave (HIFU) technology to interfere in the formation of subcutaneous fat, at the same time, to stimulate autolysis of triglycerides in lipocyte, thus, achieving slimming effects.


5 Vaginal Hifu Handle




- 360° rotation emission: total care for vaginal.
- Non-invasive, no downtime, no recovery time.
- It is easy and convenient for operation, and no consumable items are required.

Potent firming effect to improve sexual life
Private lubrication of increasing secretion and eliminating dryness
Improving private sensitivity by cell regeneration for reaction enhancement
Comprehensive improvement of private health and reducing infection

5D-HIFU Interface



Brand



Handpieces



Parameter Interface



Parameter Interface

Clients' feedback

