


360 FROZEN SLIMMING FOUR-IN-ONE BEAUTY APPARATUS

- ✓ 40K cavitation
- ✓ Body RF
- ✓ Face RF
- ✓ Lipo laser




NEW TECHNOLOGY---CRYO ANGLE

Cryo Angle is the latest fat freezing cooling technology which uses a special 360° applicator to target stubborn fat that is resistant to changes in diet and exercise, effectively freezing, destroying, and permanently eliminating the fat cells beneath the skin without damaging the surrounding layers.



CONTRASTING CASE

Develop the most excellent and innovative products, Constantly go beyond, the pursuit of perfection, Honesty, innovation as the soul



Frozen weight loss

Four handles can work at same time

- Temperature -10°C/+45°C
- 360°Frozen Fat-Dissolving Series
- Handle switch/screen synchronization

DISSOLUTION OF FAT CELL

FAT CELL WILL BE DECOMPOSED IN CRYOGENIC ENVIRONMENT, AND ELIMINATED BY METABOLISM, SO THE FAT LAYER WILL BE THINNER AND THINNER



- Apply Coolshaping applicator at the target area. Start cryo-cooling to reach a target temperature.
- The fat cells in the target area begin responding to cryo-crystallization.
- Following the treatment, fat cells reach an biological death and removed by metabolism process.
- The fat cells are eliminated.

1 High quality Working Heads

EXCHANGEABLE CRYO HANDLE CUPS

EASY TO CHANGE CRYO HANDLE CUPS

- CONTOURS ARE SECURELY FASTENED TO THE APPLICATOR
- QUICKLY CHANGE CONTOURS FOR YOUR NEXT TREATMENT

LARGE SIZE
Large cryo handle size: 1#195*85mm; 2#210*90mm; 3#215*100mm; 4#230*105mm

#3/#4 APPLICATOR
for flanks, curved cup design allows for better placement and fit on curved or narrow parts of the body.

#2 APPLICATOR
best matches with the contour of abdomen, most commonly used applicator used for this area.

MEDIUM SIZE
Small cryo handle size: 1#140*70mm; 2#155*80mm; 3#165*90mm

#1 APPLICATOR
ideal for longer vertical areas of fat (I.E. abdomen and flanks), also suggested for the arms.

COOLING CUP

Front Chin, Arm/Pit Fat, Knees

Back Chin, Arm/Pit Fat, Knees

Mini freeze
HANDLE SIZE: 80*40mm

- ★ Ideal for smaller bulges
- ★ Safety targets delicate, hard-to-reach areas
- ★ 35 minute cycle, Complete results in 1-2 visits
- ★ A smooth, contoured, youthful jawline and chin
- ★ No scalpel, injections, incisions, traumatic suction, or anesthesia

FreezeCurve

- Applicator of choice for flanks
- Curved cup design allows for better placement and fit on curved or narrow parts of the body

FreezeCore

- Best matches with the contour of abdomen
- Most commonly used applicator used for this area

FreezeFit

- Ideal for longer vertical areas of fat (I.e. abdomen and flanks)
- Also suggested for the arms and inner thighs (OUS)

Large size

Medium size

2 Product advantages

360 ANGLE SURROUNDING CRYOLIPOLYSIS SLIMMING MACHINE

40K Cavitation
direct into the fatty layer, speedily vibrate deep-seated cellulite, produce numberless vacuum cavitation, mightly strike the fatty cells, let them produce inner cracking, and dissolve to be the free fatty acid.

Body RF
uses proprietary technology to create a highly efficient and tightly woven energy matrix. This dense energy matrix penetrates multiple layers of the dermis heating it from the inside out. Complete skin penetration ensures exponentially improved clinical efficacy and products clearly visible results.

Face RF
RF radio frequency is one kind of high frequency electromagnetic wave. It can penetrate the dermis and subcutaneous tissue. Face RF can be used for face lifting, wrinkle removal.

Lipo laser
utilizes the latest low level cold laser technology to reduce girth and spot fat over almost all areas of the body without surgery, downtime or redness. A typical 40-minute waistline treatment can reduce an individual by 1/2 to 3/4 of an inch while a full protocol of eight treatments usually results in multiple inch loss.

CRYO 1, CRYO 2, CRYO 3, CRYO 4, CRYO 5, Lipo Laser, 40K, BodyRF, FaceRF

Technical Specifications:

Voltage	AC110V/220V 50-60Hz
Power	800W
Negative pressure	0-80Kpa
Screen	10.4 inches touch screen
RF	5Mhz
Cryo Quantity	5 pcs
Cryolipolysis	-10to -45°C
Cryolipolysis handle	Large/Medium/Small
Coolant	Purified water or special coolant Ambient

Dimensions: 45cm x 62cm x 170cm